



F2 Adjustable-Volume Pipetters

Cat. No.	Range	Increments	Accuracy	Precision	Compatible Tips	Color Code
4642010	0.2-2 μ L, micro	0.002 μ L	\pm 12.0 to 2.5%	10.0 to 2.0%	Flex 10, 10, 50	Pink
4642020	0.5-5 μ L, micro	0.01 μ L	\pm 6.0 to 1.5%	5.0 to 1.0%	Flex 10, 10, 50	Pink
4642030	1-10 μ L, micro	0.02 μ L	\pm 2.5 to 1.0%	2.0 to 0.5%	Flex 10, 10, 50	Pink
4642040	1-10 μ L	0.02 μ L	\pm 3.5 to 1.0%	3.0 to 0.8%	Flex 200, 250 Univ.	Yellow
4642050	2-20 μ L, micro	0.02 μ L	\pm 3.0 to 1.0%	2.5 to 0.4%	50 micro	Turquoise
4642060	2-20 μ L	0.02 μ L	\pm 3.0 to 1.0%	2.5 to 0.4%	Flex 200, 250 Univ.	Yellow
4642120	5-50 μ L, mico	0.1 μ L	\pm 0.3 to 0.15%	2.5 to 0.3%	50 micro	Turquoise
4642130	5- 50 μ L	0.1 μ L	\pm 0.3 to 0.15%	2.5 to 0.3%	Flex 200, Flex 300, 250 Univ., 200 Ext.	Yellow
4642070	10-100 μ L	0.2 μ L	\pm 3.0 to 0.8%	1.0 to 0.2%	Flex 200, 250 Univ., 300, 200 Ext.	Yellow
4642080	20-200 μ L	0.2 μ L	\pm 1.8 to 0.6%	0.7 to 0.2%	Flex 200, 250 Univ., 300, 200 Ext.	Yellow
4642090	100-1000 μ L	1 μ L	\pm 1.0 to 0.6%	0.6 to 0.2%	Flex 1000, 1000, 1000 Ext.	Blue
4642100	0.5-5 mL	0.01 mL	\pm 2.0 to 0.5%	0.8 to 0.2%	5 mL	Green
4642110	1-10 mL	0.02 mL	\pm 2.0 to 0.5%	0.8 to 0.2%	10 mL, Flex 10 mL Ext.	Red